

BOUNDER *Plus*

High Volume Low Pressure Off-Road Package



Available on
**BOUNDER *Plus* H-Frame
Power Wheelchairs**

Travels Through:

- * *Grass*
- * *Dirt*
- * *Gravel*
- * *Sand*
- * *Snow*



True Off-Road Performance

21st Century SCIENTIFIC, Inc.

4931 N Manufacturing Way, Coeur d'Alene, ID 83815
(208) 667-8800 FAX (208) 667-6600
21st@wheelchairs.com

www.wheelchairs.com
(800) 448-3680

BOUNDER *Plus*

High Volume Low Pressure Off-Road Package (HVLP OR)



Increased ground clearance under battery box

Power wheelchair users need their chairs. This is a fact. Good seating, smooth ride, and reliability are necessary. Sometimes it is also NECESSARY to drive "off-road" and perhaps your environment is not all paved city streets. When the going gets tough, TOUGH BOUNDER users really get going - going with the OFF-ROAD OPTION that makes it all possible!*



14" x 5.40" HVLV tires and off-road wheelie wheels



Chain drive

The special 14" x 5.40" low pressure pneumatic tires provide a very smooth ride. Going over rough terrain is no longer bumpy and jarring. The aggressive knobby tire design provides excellent traction when traversing loose terrain, mud or even snow. The Super Torque 4 medium speed motors have all the torque needed to climb steep hills. The Off-Road Option Package is what power wheelchair users want in a performance chair.

High Volume Low Pressure Off-Road Package (HVLP OR) adds these features to your **BOUNDER *Plus*** Power Wheelchair:

- * **14" x 5.40" Low Pressure Pneumatic Knobby Black Tires** (ATV Style) - provides ultimate grip and a smooth, comfortable ride
- * **9" x 3.50" Black Knobby Pneumatic Caster Tires** - foam filled tires available
- * **Super Torque 4 Medium Speed Motors** - provides superb torque and a 6.5 mph top speed
- * **Available with choice of chain or Panther drive**
- * **4.5" Ground Clearance Under Battery Box** - 1" extra ground clearance to allow for easier transition over obstacles
- * **Off-Road Wheelie Wheels** - provides greater protection against tipping when traveling up hills in loose soil
- * **28.5" Overall Width** - from outside of drive wheel to outside of drive wheel

Add other options to this package for a full featured chair. For example, add our stable seat elevator or a tilt system.



Bounder *Plus* H-Frame with Panther drive, PSE2 power seat elevator and 50 degree power tilt, shown in the elevated position

* Never use your power wheelchair in unsafe conditions. Users should not attempt off-road driving alone and caution is advised. Make sure you have someone to help you in case you get into trouble. Avoid unsafe hills, inclines or obstacles that may cause your chair to tip over or become stuck. Wear proper positioning straps, helmet and goggles if needed. Pack of sand and loose soil will affect performance.

21st does not condone reckless driving of your power wheelchair and will not be held liable for damage, injury or death as a result. User accepts all liability if driven unsafely or in unsafe conditions.

21st Century Scientific, Inc. reserves the right to make changes at any time without notice.

Issue date: June 2009