

New



Yes, you can.

➤ Invacare® **Action® 2NG**

The Next Generation has arrived!

The all-new **Action 2NG** has been developed from the tried-and-tested Action 3NG concept. With new customization features such as depth-adjustable seat upholstery, the specifications meet the basic needs of the users while offering greater adjustability. The next generation **Action 2NG** shares the same essential platform as the Action 3NG, offering a low total weight for easy handling combined with good driving abilities. Due to the shared use of common components, the **Action 2NG** fits perfectly into the Action family.



Light and comfortable

With a weight of just 15kg, the Invacare® **Action 2NG** fulfills all essential functional requirements. The stable folding frame means the wheelchair is easy to handle, and the light weight makes transportation trouble-free. Cross-compatibility with other members of the Action family ensures that the **Action 2NG** can be easily adapted to meet individual user needs.

Quality of life

The tried-and-tested technology of the Action family products allow the user to go wherever they desire. Cross-compatibility with Action 3NG make this an ideal, basic, fleet product.



Features



NEW Depth-adjustable seat upholstery

Decrease seat depth by 50mm for customised seat size.



Swing in and out footrest


The 80° footrest swings in and out for easy transfers in or out of the chair.



Adjustable center of gravity

The newly-designed rear wheel plate can be adjusted by 5 cm.

Frame Colour



Anthracite Black

Upholstery



Black Nylon

Please note colors may vary slightly from those displayed above.

Technical Data

41 cm, 45 cm, 48 cm	40 - 45 cm	48.5 cm	43 cm
26 cm	35 - 48 cm	80° / 0° - 80°	7° / 0 - 30°
58.5 - 68.5 cm	93.5 cm	82 / 84.5 cm	41 cm
			Warranty: 2 years frame & crossbraces. Limited 1 year components.
16.2kg Self Propel 14.9kg Transit	max. 125 kg	Crash tested ISO 7176-19*	

* Our wheelchairs comply above ISO standards and with ISO 10542, and have additional features providing increased levels of occupant security and safety whilst travelling in a motor vehicle. Wheelchair users should however transfer to the vehicle seat and use the vehicle installed restraint system whenever it is feasible.



Yes, you can.

Your Local Dealer:

All information quoted is believed to be correct at time of print. Invacare® reserves the right to alter product specifications without prior consultation.

Invacare New Zealand
4 Westfield Place (PO Box 62-124) Mt Wellington Auckland, New Zealand
Freephone: 0800 INVACARE (0800 468 222) Freefax: 0800 807 788
Email: sales@invacare.co.nz Web: www.invacare.co.nz

Invacare Australia Pty Ltd
1 Lenton Place (PO Box 5002) North Rocks NSW 2151, Australia
Freephone: 1800 460 460 Fax: 1800 814 367
Email: sales@invacare.com.au Web: www.invacare.com.au

Invacare Asia Ltd
1 Lenton Place (PO Box 5002) North Rocks NSW 2151, Australia
Phone: +61 2 8839 5330 Fax: +61 2 8839 5331
Email: asiassales@invacare.com Web: www.invacare.com